

ambulate

CLINICAL TRIAL

VASCULAR CLOSURE FOR EP PROCEDURES

Proven by EPs

VASCADE MVP was evaluated in the AMBULATE pivotal trial. Natale et al. AMBULATE was a multicenter, randomized clinical trial with 204 patients randomized 1:1 to VASCADE MVP or manual compression

VASCADE MVP®
VENOUS VASCULAR CLOSURE SYSTEM



vs.



**MANUAL
COMPRESSION**



204
patients



SEP
2018

OCT

NOV

DEC

JAN

FEB

MAR
2019



13
sites



28
investigators

The only
FDA-approved
closure device
for use
following cardiac
ablations.

Time to Ambulation



Reduction in Median
Time to Ambulation¹

Patient Satisfaction



Improvement in
Patient Satisfaction¹

Opioid Use



Reduction in Opioid Use¹

Designed for EP Procedures

- Simple and Easy to Use
- Enables Multiple Site Closures in Close Proximity
- Nothing Permanent Left Behind for Future Procedures

VASCADE MVP®

VENOUS VASCULAR CLOSURE SYSTEM

**Early Ambulation.
Simple. Proven.**

NOW FDA APPROVED



Indicated for 6-12F
Inner Diameter Venous
Vascular Closure

The Only Closure Device Indicated for Multi-site, Multi-limb Femoral Venous Closure



AMBULATE Pivotal Trial – Summary of Endpoint Analysis

Attribute		VASCADE MVP N=199	Manual Compression N=209	p-value	
Safety Endpoint	Major Complications	0%	0%	–	
	Minor Complications	1.0%	2.4%	0.45 ⁴	
		N=100	N=104		
Efficacy Endpoint*	Time to Ambulation (TTA)	Mean	2.8 hours	6.1 hours	<0.0001 ²
		Median	2.2 hours	6.1 hours	<0.0001 ²
	Time to Closure Eligibility (TTCE)	10.5 mins	37.6 mins	<0.0001 ²	
	Total Post-Procedure Time (TPPT)	3.1 hours	6.8 hours	<0.0001 ²	
	Time to Discharge Eligibility (TTDE)	3.1 hours	6.5 hours	<0.0001 ²	
	Time to Discharge (TTD)	21.8 hours	21.8 hours	0.98 ²	
	Procedure Success	98%	99%	–	
			N=369	N=382	
	Time to Hemostasis (TTH)	6.1 mins	13.7 mins	<0.0001 ²	
			N=363	–	
Device Success	97%	N/A	–		

AMBULATE Pivotal – Patient Satisfaction Scores⁵

Attribute	VASCADE MVP N=100	Manual Compression N=102	Nominal p-value ³	Outcome	
Patient-Reported Satisfaction Scores	All Patients				
	Duration, current episode	8.3±2.4	5.1±3.4	<0.0001	Scores 63% higher for MVP
	Discomfort, current episode	7.2±3.1	5.3±3.1	<0.0001	Scores 36% higher for MVP
	Pain, current episode	7.5±3.2	6.0±3.4	0.001	Scores 25% higher for MVP
	Patients with a previous cardiac ablation procedure				
	Duration, compared to previous	7.9±2.3	5.6±3.0	0.001	Scores 41% higher for MVP
	Discomfort, compared to previous	7.5±2.1	5.4±2.8	0.001	Scores 39% higher for MVP
Pain, compared to previous	7.7±2.8	5.5±2.9 (N=38)	0.002	Scores 40% higher for MVP	

AMBULATE Pivotal Results – Pain Medications & Opioid Use

Pain Medication Usage	VVCS (N=100)		Manual Compression (N=104)		Nominal p-value ⁴
Yes	24	24%	51	49%	0.0003
Reduction in Overall Pain Meds			51%		
Opioid Usage	VVCS (N=100)		Manual Compression (N=104)		Nominal p-value ⁴
Yes	15	15%	37	36%	0.001
Reduction in Opioid Use			58%		

1. AMBULATE Clinical Study Report

2. p-values from two-sided t-test for means, and two-sided Wilcoxon rank sum test associated with medians, unadjusted for stratification factor

3. Nominal p-values by two-sided t-test, not pre-specified, however actual results imply superiority

4. Nominal p-value by 2-sided Fisher's exact test

5. Patient satisfaction surveys administered prior to discharge. Rated on scale of 1-10, with "10" being most satisfied. Data are presented as mean ± SD.

* Mean timepoints shown, unless otherwise noted



P: 866.602.6099 F: 866.602.1795
customerservice@cardivamedical.com
www.CardivaMedical.com